

Stainless Steel Dual Head Stethoscope

LAB-7100
LAB-7300

- Stainless steel dual head chestpiece with non-chill ring for patient comfort
- Internal spring binaurals ensure acoustic sensitivity
- Contoured ear tips
- 22 inch (55.88 cm) Latex Free Tubing

LAB-7100 - Adult ■

LAB-7300 - Infant ■
LAB-7300 shown



Contact your medical supply dealer for the complete line of Graham-Field diagnostic equipment.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com

DUAL HEAD STETHOSCOPES

Superior Performance

Manufactured for Durability

Unsurpassed Quality



grahamfield

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company



© 2006, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.

Color variations may occur due to monitor or printer settings.

STETH-INS-LAB-RevC24

Dual Head Stethoscope 400

- Flexible PVC tubing in a choice of colors
- Non-chill diaphragm and recessed bell retaining rim for patient comfort
- Adjustable metal binaural
- 22 inch (55.88 cm) Latex Free Tubing

400 ■
400BL ■

400GR ■
400GY ■

400P ■
400R ■
shown in pink [e]



Panascope™ Dual Head Stethoscope

500
510

- Non-chill bell and diaphragm retaining rim
- Color coordinated retaining rim, eartips, and flexible PVC tubing
- Lightweight aluminum binaural
- Includes additional diaphragm and set of eartips
- 22 inch (55.88 cm) Latex Free Tubing

500 - Adult ■ shown in gray [e]
500GY - Adult ■



Sprague Rappaport Stethoscope 602

602

- Five interchangeable chestpieces (two diaphragms, 3 bells) facilitates detection of high, medium and low frequency sounds
- Ideal for infant, pediatric, and adult use
- Includes 3 eartip sets to maximize comfort
- 22 inch (55.88 cm) Latex Free Tubing

602 ■
602BL ■

602R ■
602TL ■
shown in teal [t]

