



BATTERY POWERED LIFTS

PRO BATTERY-POWERED FLOOR LIFT

Lumex floor lifts are designed to provide safer and more secure transfers of patients by caregivers. Lumex lifts are very sturdy and easy to use.

Features

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Swan neck enables outstanding lifting range of 55.5" (140.97 cm)
- Six point spreader bar with 360° rotation accommodates 2-, 4-, or 6-point slings
- Large, easy to grip padded handles for easy maneuverability
- Powered high/low and base open/close simplifies operation
- Ergonomic hand pendant with hook
- Premium electronics with LCD display. Lift can be operated via the control box or the hand pendant
- Emergency stop button and manual lowering mechanism
- Alarm and visual indicator on control box and hand pendant for low battery charge (24V DC motor) or overload
- Battery may be charged while on the lift or remotely in charging caddy
- Spare battery and charging caddy included
- Durable twin wheel casters, 3" (7.62 cm) front, 4" (10.16 cm) rear braking
- Designed to meet the requirements of HCPCS Code E0635
- How-to-use video included



Warranty

- Frame and spreader bar – Three (3) year limited warranty
- Actuator – Two (2) year limited warranty
- Batteries – One (1) year limited warranty

	Floor Clearance		Spreader Bar Height (measured from sling attachment)				Base Width (outer measurement)				overall Length		product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg			
LF500	4.3	10.92	19.3	49.02	74.8	189.99	25.6	65.02	39.4	100.08	49.2	124.97	116.8	53.	500	226.8	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-1600038Rev1B23