



BATTERY POWERED LIFTS

BATTERY-POWERED FLOOR LIFT


Lumex floor lifts are designed to provide safer and more secure transfers of patients by caregivers. Lumex lifts are very sturdy and easy to use.

Features

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Six point spreader bar with 360° rotation accommodates 2-, 4-, or 6-point slings
- Large, easy to grip handles for easy maneuverability
- Intuitive foot pedal leg opening mechanism frees up hands
- Ergonomic hand pendant with hook
- Emergency stop button and manual lowering mechanism
- Warning beep and visual indicator on control box for low battery charge (24V DC motor)
- Accessory storage bag and how-to-use video included
- Optional digital scale DSC250
- Designed to meet the requirements of HCPCS Code E0635

Warranty

Frame and spreader bar – Three (3) year limited warranty
 Actuator – Two (2) year limited warranty
 Batteries – One (1) year limited warranty



	Floor Clearance		Spreader Bar Height (measured from sling attachment)				Base Width (outer measurement)				overall Length		product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	minimum inch	maximum cm	minimum inch	maximum cm	minimum inch	maximum cm	inch	cm	lb	Kg	lb	Kg			
LF1050	4.7	11.94	27.5	69.85	68.1	172.97	40.2	102.11	26	66.04	37.8	96.01	114.5	51.9	400	181.4	1 ea



GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
 Lumex is a registered trademark of GF Health Products, Inc.