



BATTERY POWERED LIFTS

EASY LIFT - BARIATRIC PATIENT LIFTING SYSTEM

Lumex floor lifts are designed to provide safer and more secure transfers of patients by caregivers. Lumex lifts are very sturdy and easy to use.

Features

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Six point spreader bar with 360° rotation accommodates 2-, 4-, or 6-point slings
- Unique, ergonomic foot pedal leg opening
- Large, easy to grip handles for easy maneuverability
- Lower base height maximizes stability and fits under most beds
- Rubber base coating protects furniture and walls
- 24V DC motor allows quiet smooth operation
- Ergonomic hand pendant with hook
- Warning beep for low battery charge and weight capacity levels
- Easy access emergency stop button
- Emergency manual lowering device
- White powder-coated finish
- Accessory storage bag included for slings and charger
- Optional digital scale available
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS Code E0635



Warranty

Frame and spreader bar – Three (3) year limited warranty
 Actuator – Two (2) year limited warranty
 Batteries – One (1) year limited warranty



	Floor Clearance		Spreader Bar Height (measured from sling attachment)				Base Width (outer measurement)				overall Length		product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	minimum	maximum	minimum	maximum	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
LF1090	4.7	11.94	27.2	69.09	70.9	180.09	50	127.	26	66.04	43	109.22	137	62.1	600	272.2	1 ea



GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
 Lumex is a registered trademark of GF Health Products, Inc.