



## BATTERY POWERED LIFTS

# BATTERY-POWERED SIT-TO-STAND LIFT

Sit-to-Stand lifts are designed for patients who have difficulty standing up on their own from a seated position.

Patients must be able to actively participate in the lifting process and be able to partially support their weight.

### Features

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Spreader bar provides various hand-hold locations for the patient and accommodates various sling styles
- Large, easy to grip handles for easy maneuverability
- Intuitive foot pedal leg opening mechanism frees up hands
- Adjustable knee pad unit
- Quick-release removable foot platform for gait training
- Ergonomic hand pendant with hook
- Emergency stop button
- Warning beep and visual indicator on control box for low battery charge (24V DC motor)
- Accessory storage bag and how-to-use video included
- Designed to meet the requirements of HCPCS Code E0635

### Warranty

Frame and spreader bar – Three (3) year limited warranty  
 Actuator – Two (2) year limited warranty  
 Batteries – One (1) year limited warranty



	Floor Clearance		Spreader Bar Height (measured from sling attachment)				Base Width (outer measurement)				overall Length		product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	minimum	maximum	minimum	maximum	minimum	maximum	inch	cm	inch	cm	lb	Kg	lb	Kg	
LF2020	4.7	11.94	36.4	92.46	70.8	179.83	25.5	64.77	37.8	96.01	40.4	102.62	118	53.5	400	181.4	1 ea



GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
 The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
 Lumex is a registered trademark of GF Health Products, Inc.