



## BATTERY POWERED LIFTS

# BARIATRIC EASY LIFT STS

### Features

- Heavy-gauge steel construction
- Designed to provide quick, safe assistance to those who have trouble standing
- Unique, ergonomic foot pedal base opening
- Large, easy to grip handles for easy maneuvering
- Lower base height maximizes stability and fits under most beds
- Rubber base coating protects furniture and walls
- Adjustable knee pad unit
- Quick release removable foot platform
- Optional buttock strap for patients requiring additional support
- 24V DC motor allows quiet, smooth operation
- Ergonomically designed hand pendant with pendant hook
- Warning beep for low battery charge and weight capacity levels
- Easy access emergency stop button
- White powder-coated finish
- Accessory storage bag included for slings and charger
- Interchangeable battery works with other Lumex Easy Lifts
- Easy access emergency stop button
- Maximum Weight Capacity  
600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS Code E0635

### Warranty

- Frame and spreader bar – Three (3) year limited warranty
- Actuator – Two (2) year limited warranty
- Batteries – One (1) year limited warranty



	Floor Clearance		Spreader Bar Height (measured from sling attachment)				Base Width (outer measurement)				overall Length		product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
			minimum	maximum	minimum	maximum	lb	Kg	lb	Kg							
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg			
LF2090	4.7	11.94	36.4	92.46	70.8	179.83	25.5	64.77	37.8	96.01	40.4	102.62	134	60.8	600	272.2	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF1200084RevB23