

## BATTERY POWERED LIFTS

# POWERED SIT-TO-STAND 500

Sit-to-Stand lifts are designed for patients who have difficulty standing up on their own from a seated position. Patients must be able to actively participate in the lifting process and be able to partially support their weight

### Features

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Sling hooks on spreader bar facilitate easy placement of sling shoulder straps
- Premium electronics with LCD display. Lift can be operated via the control box or the hand pendant
- Large, easy to grip handles for easy maneuverability using 4" locking and 3" non-locking casters
- Powered high/low and base open/close simplifies operation (STS500 only)
- Intuitive foot pedal leg opening mechanism frees up hands for manual base option (STS500M only)
- Adjustable knee pad unit
- Quick-release removable foot platform for gait training
- Ergonomic hand pendant with hook
- Emergency lowering device on actuator shaft can be used in case of power failure
- Warning beep and visual indicator on control box for low battery charge (24V DC motor)
- Maximum Weight Capacity of 500 lb (226.8 kg) EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS Code E0635

### Warranty

**Frame and spreader bar:** Three (3) year limited

**Actuator:** Two (2) year limited

**Batteries:** One (1) year limited



Specification		STS500	STS500M
External base width	Open	38.6" (98 cm)	39.4" (100 cm)
	Closed	26" (66 cm)	25.6" (65 cm)
Internal width, lifting boom		18.8" (47.8 cm)	
Base height (for clearance)		4.3" (11 cm)	
Base legs length		40.6" (103 cm)	
Overall base length		40.6" (103 cm)	
Lifting height range as measured from sling hooks	Minimum	36.2" (92 cm)	19.3" (49 cm)
	Maximum	70" (178 cm)	74.8" (190 cm)
	Range	33.8" (86 cm)	55.5" (141 cm)
Mast height from floor		44.1" (112 cm)	
Pad height from floor		13.8" (35 cm)	
Maximum weight capacity		500 lb (230 kg), EVENLY DISTRIBUTED	
Caster size		4" diameter locking; 3" diameter non-locking	
Lifting mechanism		24 DC actuator	Manual
Battery charger		AC 100V-240V to DC 29V, 1.5A	N/A



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2023, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.