

EXAM ROOM SUPPLIES

EMESIS BASINS

Features

- 18-8 (304) stainless steel with seamless construction.
- Shallow basin of curved, kidney-shaped design
- Autoclavable



		overall Length		overall Depth		overall Height		Volume		UOM
		inch	cm	inch	cm	inch	cm	oz	ml	
3231	Small	8.125	20.6	3.625	9.2	1.75	4.4	13	384.5	1 ea
3232	Large	10	25.4	4.5	11.4	2.125	5.4	26	768.9	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.