

EXAM ROOM SUPPLIES

DRESSING JARS WITH LIDS

Features

- Recessed lid handles
- Thickness: 0.7mm
- 18-8 (304) stainless steel with seamless construction



	overall Height		overall Width		Volume		UOM
	inch	cm	inch	cm	qt	l	
3233	6.875	17.5	4.875	12.4	2.125	2	1 ea
3233-1	5	12.7	4	10.2	1	0.9	1 ea
3233-2	4.125	10.5	2.5	6.4	0.5	0.5	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.