

EXAM ROOM SUPPLIES

STAINLESS STEEL SOLUTION BOWL

Features

- Rounded lip for easy pouring
- Capacity: 7 qts
- Seamless stainless steel
- Autoclavable



| | Height | | Diameter | | Volume | | UOM |
|------|--------|------|----------|------|--------|-----|------|
| | inch | cm | inch | cm | qt | L | |
| 3247 | 4.625 | 11.7 | 13.625 | 34.6 | 7 | 6.6 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.