

EXAM ROOM SUPPLIES

# STAINLESS STEEL INSTRUMENT TRAYS

### Features

- Designed to hold medical supplies
- Flat bottom tray
- 18-8 (304) stainless steel
- Sterilizable



Optional lid with recessed grip



Optional lid with easy to grip raised strap handles

		overall Length		overall Width		overall Height		UOM
		inch	cm	inch	cm	inch	cm	
3254	With Strap Handle Cover	8.5	21.6	3	7.6	1.5	3.8	1 ea
3255	With Strap Handle Cover	8.875	22.5	5	12.7	2	5.1	1 ea
3256	With Recessed Grip Cover	12.125	30.8	7.625	19.4	2.125	5.4	1 ea
3257	With Strap Handle Cover	17	43.2	4	10.2	2.125	5.4	1 ea
3258	Without Cover	8.875	22.5	5	12.7	2	5.1	1 ea
3259	Without Cover	12.25	31.1	7.625	19.4	2.125	5.4	1 ea



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved. Graham-Field is a registered trademark of GF Health Products, Inc.