

LAMPS

B-SERIES BATTERY-OPERATED UV HAND-HELD LAMP

Features

- Anodized aluminum housing
- Aluminum reflector optimizes UV irradiance for brighter fluorescent response
- Suitable for long wave lamp application
- Highest UV irradiances of any units of its type
- Compact and lightweight: requires 5 C-batteries (not included)
- Can be operated for approximately five hours on fresh batteries
- NiCad rechargeable batteries can be used
- Transparent plastic tube guard to ensure against accidental breakage
- Tube: 4-Watt (BLE-22B)
- Typical Peak UV Intensity at six inches: 730 $\mu\text{W}/\text{cm}^2$
- Made in USA, except bulbs

USED TO DIAGNOSE only



Tube	Typical Peak UV Intensity at 6" (15.2 cm)		overall Length		overall Width		overall Height		UOM
	$\mu\text{W}/\text{cm}^2$		inch	mm	inch	mm	inch	mm	
2218	4-Watt (BLE-22B)	750	8	20.3	3.25	8.3	2.5	6.4	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.