

EXAM ROOM SUPPLIES

FLAT TYPE INSTRUMENT TRAYS

Features

- Raised edge all around the perimeter
- Accommodates supplies, instruments, and tools
- 18-8 (304) stainless steel with seamless construction



| | overall Length | | overall Width | | overall Height | | UOM |
|------|----------------|------|---------------|------|----------------|-----|------|
| | inch | cm | inch | cm | inch | cm | |
| 3261 | 13.625 | 34.6 | 9.75 | 24.8 | 0.625 | 1.6 | 1 ea |
| 3262 | 15.125 | 38.4 | 10.5 | 26.7 | 0.625 | 1.6 | 1 ea |
| 3264 | 19 | 48.3 | 12.5 | 31.8 | 0.625 | 1.6 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.