

EXAM ROOM SUPPLIES

# STAINLESS STEEL KICK BUCKET/STAND SET

**Features**

- Bucket pours easily without removal from stand
- Vinyl covered ring on stand prevents marring of walls and floors
- Smooth rolling 2" ball bearing conductive casters
- Capacity: 12 ½ quarts
- Bucket dimensions: 11 5/8" D x 9 1/4" H
- Stainless steel



		bucket Height		bucket Diameter		bucket Volume		caster Diameter		UOM
		inch	cm	inch	cm	qt	L	inch	cm	
3267	Kick Bucket and Stand Set	9.25	23.5	11.625	29.5	12.5	11.8	2	5.1	1 ea
3267P	Utility Pail Only	9.25	23.5	11.625	29.5	12.5	11.8	—	—	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2011, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.