

EXAM ROOM SUPPLIES

# STAINLESS STEEL DRY HEAT STERILIZER

The improved technique for dry heat sterilization in the shortest possible time. Sterilizes instruments, needles, syringes, dressings at 320 – 360° F (160 – 180° C). Sharp instruments remain sharp; no steam or moisture to corrode the finest edge. Exclusive internal-wall tubular heating system is burn-out proof. 115V AC. UL/CSA listed.

**The Best Method for Sharp Instruments**

Dry Heat Sterilizer kills bacteria without moisture and will not dull sharp instruments. Plus, they can be used to dry wet packs or as a back-up sterilizer in your practice - the perfect complement to an autoclave.

**Maintenance-Free**

Unlike wet heat methods, Dry Heat Sterilizers require no routine cleaning.

**Safe and Simple to Operate**

Dry heat sterilization works without pressure, steam, or chemicals, which increases safety.

**Convenient**

Utility trays can be sterilized with lids on and stored without contamination. Also, they have automatic timers with heat "hold" capability for continuous use in a busy practice.

**High-Quality**

Stainless steel construction means no rust.



		UOM
8375	500 Watt Dry Heat Sterilizer with 2 Trays	1 ea
8375-220	500 Watt Dry Heat Sterilizer with 2 Trays (220V)	1 ea
8376	735 Watt Dry Heat Sterilizer with 3 Trays (Hospital Model)	1 ea




**Dry Heat Self-Seal Bags**

- Made of medical grade paper with an indicator strip
- Brown when exposed to high temperature  
2.5' x 1.5' x 10.5" (6.4 cm x 3.8 cm x 26.7 cm)

UOM		
8375-250	Dry Heat Self-Seal Bags	250 ea/bx

SPORE TEST STRIPS: The Centers for Disease Control (CDC) recommend the sterilizer be tested weekly. For Dry Heat Sterilizers, use Bacillus Atrophaeus spore test strips.\*  
\*Not included

GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700  
Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.