

EXAM ROOM SUPPLIES

ULTRASOUND GEL WARMER

Features

- White, metal gel warmer holds one 8 oz (236.6 ml) bottle
- The new and improved sleeve holds the bottle deeper than the old style allowing for more even heating of the gel
- Heats gel from 90°F- 120°F (+/- 3°F) (32.2°C – 48.9°C)



		Temperature Range		UOM
		°F	°C	
GF108	1 bottle capacity	90 – 120	32.2 – 48.9	1 ea
GW208	2 bottle capacity			1 ea

-GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2005, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.