


DIAGNOSTIC INSTRUMENTS

# CRILE HEMOSTATIC FORCEPS

## Features

- Serrations along the entire length of the jaw
- Ratchets used to lock the forceps and hold blood vessels
- Stainless steel
- Professional use only



		Length		UOM
		inch	cm	
2676	Straight	5.5	14	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.