

DIAGNOSTIC INSTRUMENTS

MAYO HEGAR NEEDLE HOLDER

Features

- Ratchets lock the needle holder with controlled grip and smooth grip release
- Stainless steel
- Professional use only



	Length		UOM
	inch	cm	
2708	5.5	14	1 ea
2709	6	15.2	1 ea
2710	7	17.8	1 ea
2711	8	20.3	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.