

DIAGNOSTIC INSTRUMENTS

# HALSEY NEEDLE HOLDER

## Features

- Ratchets used to lock the needle holder with controlled grip and smooth grip release
- Smooth Jaw
- Stainless steel
- Professional use only



|      | Length |      | UOM  |
|------|--------|------|------|
|      | inch   | cm   |      |
| 2713 | 5      | 12.7 | 1 ea |



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.