

DIAGNOSTIC INSTRUMENTS

# BACKHAUS TOWEL FORCEPS

## Features

- Used for securing towels and drapes
- Sharp overlapping pointed tips
- Ratchets used to lock the forceps
- Reusable
- Medical-grade and premium quality forceps
- Stainless steel



	Length		UOM
	inch	cm	
2725	3.5	8.9	1 ea
2726	5.25	13.3	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.