

DIAGNOSTIC INSTRUMENTS

# FOERSTER SPONGE FORCEPS

## Features

- Used to hold or mobilizing sponge to absorb excess fluids and blood
- Oval tips with serrations
- Ratchet handle helps to hold the sponges firmly
- Reusable
- Stainless steel



		Length		UOM
		inch	cm	
2732	Smooth	9.5	24.1	1 ea
2733	Serrated	9.5	24.1	1 ea
2734	Curved, Serrated	9.5	24.1	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.