


DIAGNOSTIC INSTRUMENTS


# POCKET GONIOMETER

## Features

- 6" (15.2 cm) clear plastic goniometer
- The orthopedic goniometer is used to measure a joint's range of motion (ROM). There are two "arms"—one that is stationary and one that is movable—that are hinged together
- One arm has centimeter and inch scale



13630



13635

		overall Length		UOM
		inch	cm	
13630	360° protractor head	6	15.2	1 ea
13635	180° protractor head	6	15.2	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2008, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.