



EXAM ROOM SUPPLIES

JAEGER READING EYE CHART

The Jaeger eye chart is used for reading up close and for determining your near vision. When reviewing the chart, you will see the notation 1 next to the paragraph with the smallest text and each progressive paragraph of larger text is noted with an increase in the number. As you progress to larger lettered paragraphs, the lettering size increases for lesser visual acuity.

Features

- Printed on both sides
- Designed to be read from a distance measured in meters

	overall Width		overall Height		test Distance		UOM
	inch	cm	inch	cm	inch	cm	
1242	4.75	12.1	7	17.8	14	35.6	1 ea

graham
field

No. 7.
1.50M

able treaty, the restitution of the standards and prisoners which had been taken in the defeat of Crassus. His generals, in the early part of his reign, attempted the reduction of Ethiopia and Arabia Felix. They marched near a thou-

No. 8.
1.75M

sand miles to the south of the tropic; but the heat of the climate soon repelled the invaders, and protected the unwarlike natives of those sequestered regions

No. 9.
2.00M

The northern countries of Europe scarcely deserved the expense and labor of conquest. The forests and morasses of Germany were

No. 10.
2.25M

filled with a hardy race of barbarians who despised life when it was separated from freedom; and though, on the first

No. 11.
2.50M

attack, they seemed to yield to the weight of the Roman power, they soon, by a signal

GRAHAM-FIELD



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2005, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.