

EXAM ROOM SUPPLIES

# ROSENBAUM POCKET VISION SCREENER CARD

### Features

- Pocket size plastic eye chart
- Precise, clear, legible
- Provides 20/800 distance equivalent testing at 14"
- Visual acuity (Va) can be assessed by using the Pocket Vision Screener
- The Pocket Vision Screener gives a distance equivalent Va even though the card is held 14" from the patient
- The back of the card features a millimeter scale and an inch scale
- Plastic
- Latex free

	overall Width		overall Height		test Distance		UOM
	inch	cm	inch	cm	inch	cm	
1243-1	3.5	8.9	6.375	16.2	14	35.6	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2005, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.