

EXAM ROOM SUPPLIES

# PLASTIC VENTILATED EYE SHIELD

### Features

- Ventilated clear acrylic shield with ¼" (0.6 cm) openings for tape attachments
- Ideal for use after surgery
- Accommodates either eye
- Made in U.S.A.



		overall Dimensions		Thickness	Weight	UOM
		in	cm	mm	g	
1276-1	Plastic Ventilated Eye Shield, Adult	2.5 x 3.08	6.4 x 7.8	2	6.53	1 ea
1276-1B	Plastic Ventilated Eye Shield, Adult	2.5 x 3.08	6.4 x 7.8	2	6.53	100/bg



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.