

EXAM ROOM SUPPLIES

STANDARD FINGER RING CUTTER

Features

- Specially designed lever slips under ring to permit rotating thumb screw to cut ring
- Knurled straight handle with 7/8" rotary jeweler's saw
- Handle with care
- Professional use only
- Stainless steel
- Latex free



		Length		UOM
		inch	cm	
1372-3	Finger Ring Cutter	5.75	14.6	1 ea
1372-4	Replacement Blade	0.875	2.2	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2005, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.