

## EXAM ROOM SUPPLIES


# "ALCOHOL" IMPRINTED DISPENSER

### Features

- Polypropylene alcohol dispenser
- Capacity: 8 oz (236.6 ml)
- Features one hand dispensing
- Priming with alcohol is required prior to using the dispenser and periodically prime the dispenser throughout lifetime of the unit
- Swing lid
- Latex free



		Volume		UOM
		fl oz	ml	
1533-1	Dispenser with "Alcohol" imprinted	8	236.6	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.