

EXAM ROOM SUPPLIES

MAYO-TYPE INSTRUMENT TABLE

Features

- Height adjustable
- Stainless steel tray
- Chrome-plated steel base
- Adjustable height
- Non locking casters
- Casters increase maneuverability
- Knockdown
- Maximum Weight Capacity
25 lb - EVENLY DISTRIBUTED



| | Tray | | | | Adjustable Height | | | | Weight Capacity | | UOM |
|--------|-------|------|-------|------|-------------------|------|---------|-----|--------------------|-------|------|
| | width | | depth | | minimum | | maximum | | EVENLY DISTRIBUTED | | |
| | inch | cm | inch | cm | inch | cm | inch | cm | lb | kg | |
| GF1843 | 12.5 | 31.8 | 19 | 48.3 | 31 | 78.7 | 50 | 127 | 25 | 11.34 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.