


MEDICAL/SURGICAL | PATIENT ROOM

NARCOTIC SAFE


Features

- Intended use: safe/storage cabinet with lock
- Double door with individual key for each door
- Keys cannot be removed when lock is in open position
- Doors suspended on full length piano hinges
- High quality steel – welded construction
- Safe can be set into the wall or wall mounted
- Two adjustable removable steel shelves



	overall Height		overall Width		overall Depth		product Weight		UOM
	inch	mm	inch	mm	inch	mm	lb	kg	
3007	14.5	36.8	9.5	24.1	4	10.2	11	4.99	1 ea
3008	24	61.	16	40.6	8	20.3	28.5	12.93	1 ea

Replacement Lock & Key Sets



- Each set includes 1 lock & 2 keys

		UOM
3011	Inner Door	1 set
3012	Outer Door	1 set



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.