



BARIATRIC PATIENT SEATING

CLINICAL CARE RECLINER

Designed for patients weighing up to 850 lb (385.5 kg) when evenly distributed, the Gendron Bariatric Clinical Care Recliner features swivel and locking casters for maneuverability. High density foam Back features a 42 degree recline and Elevating Legrest support. The recline mechanism is located at the rear of the Recliner for easy caregiver access.

Product Info



Available Colors



AREAS OF USE
PACU
Rehab
Dialysis
Bariatrics
Oncology
IV Therapy
Exam Room
Patient Room
Emergency Room
Physical Therapy
Pain Management
Nocturnal Dialysis
Post Operative Care
Same Day Surgery / ASC



Specifications

Item	Description	Overall Width		Overall Depth		Overall Height		Seat Width		Seat Depth		Seat Height from floor		Back Height		Shipping ⁸ Weight		Weight ¹ Capacity EVENLY DISTRIBUTED		Seat Tilt degrees	Backrest Tilt degrees	Made in USA	Warranty in months
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg [†]	lb	kg [†]						
7150Pxxxx	Large Seat	38.25	97.16	68.00	172.72	47.50	120.65	29.25	74.30	19.25	48.90	20.00	50.80	33.00	83.82	233	105.69	850	385.55	20	42	■	36
7155Pxxxx	Medium Seat	34.25	87.00					25.25	64.14							230	104.33			20	42	■	36

[†] xxxx in the product number is the placeholder for fabric choice [†] standard rounding method used
¹ maximum weight capacity - evenly distributed - refer to user manual ⁸ approximate measurement



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.