



BATHROOM SAFETY

BARIATRIC COMMUNE

**Features**

- 750 lb (340.19 kg) weight capacity <sup>EVENLY DISTRIBUTED</sup><sup>1</sup>
- Powder-coated steel frame
- Adjustable, removable armrest
- Multiple Seat width options
- Pail is removable from the back

Product Info



**Specifications**

Item	Overall Width		Overall Depth		Overall Height Range				Seat Width		Seat Depth		Seat Height Range				Back Height above seat		Shipping Weight <sup>8</sup>				Weight <sup>1</sup> Capacity EVENLY DISTRIBUTED	
	inch	cm	inch	cm	minimum	maximum	inch	cm	inch	cm	inch	cm	minimum	maximum	inch	cm	without pallet	with pallet	lb	kg†	lb	kg†	lb	kg†
523326	31.40	79.76	22.00	55.88	19.25	48.90	23.25	59.06	26.00	66.04	20.00	50.80	20.00	50.80	24.00	60.96	12.50	31.75	50	22.68	90	40.82	750	340.19
523330									30.00	76.20							55	24.95	95	43.09				

† standard rounding method used    <sup>1</sup> maximum weight capacity - evenly distributed - refer to user manual    <sup>8</sup> approximate measurement



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved.  
Gendron is a registered trademark of GF Health Products, Inc.