

## DIAGNOSTIC INSTRUMENTS

# ADJUSTABLE FREQUENCY TUNING FORK

### Features

- Fork has adjustable masses that can be moved up and down to change the pitch
- All aluminum tuning fork, with adjustable weight scale in seven frequencies: C128, D144, E160, F170.7, G192, A213.3 and B240 (weights removed)
- Scale letters are stamped on stems



	tine Width		overall Length		UOM
	inch	cm	inch	cm	
1325	0.375	0.95	8	20.3	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.