



MEDICAL SEATING

PHYSICIAN EXAM STOOL

9000 SERIES

Lumex® offers a variety of medical stools to meet the needs of an active facility. Manufactured with the highest quality materials to ensure a long life, Lumex stools offer the most wanted features in an exam stool.

Features:

- Seat diameter comfortably fits to most users
- Pneumatic assist height adjustment for ease of use
- Choice of 3 base styles to complement facility decor
- Choice of height adjustment lever options
- Wide diameter base for greater stability
- Backrest and foot ring options
- Hard wheel casters
- Available in 24 industry standard colors to complement most decors

Product Info 9501 9702 9903 Fabric Choices



9501 Options:
black nylon composite base †

Standard Height
with Backrest
with Footring
Lab Height

9501-xx	■			
9501B-xx	■	■		
9501-FR-xx	■		■	
9501-BFR-xx	■	■	■	
9501-10BFR-xx		■	■	■

★ options available at additional cost



9702 Options:
brushed aluminum base and toe caps †

Standard Height
with Backrest
with Footring
Lab Height

9702-xx	■			
9702B-xx	■	■		
9702-FR-xx	■		■	
9702-BFR-xx	■	■	■	
9702-10BFR-xx		■	■	■

★ options available at additional cost



9903 Options:
polished aluminum base †

Standard Height
with Footring
Lab Height

9903-xx	■			
9903-FR-xx	■	■		
9903-10FR-xx			■	■

★ options available at additional cost

Specifications:

Series	Seat Diameter		Base Diameter with casters		Standard Seat Height		Lab Seat Height		Overall Height w/ Backrest				Footring Diameter		Caster Diameter hard wheel		Shipping Weight ⁸		Weight Capacity ¹ EVENLY DISTRIBUTED		Pneumatic Height Adjustment Method	Single Lever D-Handle	360° Chrome Handle release					
	inch	cm	inch	cm	low	high	low	high	low	high	low	high	inch	cm	inch	cm	lb	kg ¹	lb	kg ¹								
9501					18.50	46.99	24.00	60.96	24.50	62.23	35.00	88.90	33.00	83.82	39.00	99.06			20	9.07								
9702	16.00	40.64	24.00	60.96	17.50	44.45	23.00	58.42	23.50	59.69	34.00	86.36	32.50	82.55	37.50	95.25	18.00	45.72	2.00	5.08	22	9.98	350	158.76				
9903					18.00 ²	45.72 ²	23.00 ²	58.42 ²	22.00	55.88	32.00	81.28	18.00 ²	45.72 ²	23.00 ²	58.42 ²			23	10.43								

★ options available at additional cost † -xx in the product number is the placeholder for fabric choice † standard rounding method used
¹ maximum weight capacity - evenly distributed - refer to user manual ² overall height without backrest ⁸ approximate measurement



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved. Lumex is a registered trademark of GF Health Products, Inc.