



## MEDICAL SEATING



# PHYSICIAN SURGICAL STOOL

## 7400 SERIES

Lumex<sup>®</sup> offers a variety of medical stools to meet the needs of an active facility. Manufactured with the highest quality materials to ensure a long life, Lumex stools offer the most wanted features in an exam stool.

### Features:

- Seat diameter comfortably fits to most users
- Hydraulic height adjustment for ease of use
- Foot Control adjustment release
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate
- Wide diameter base for greater stability
- Seat rotation lock
- Silver powder coat column
- Hard wheel casters
- Crescent back and backrest options
- Available in 24 industry standard colors to complement most decors

7408 Product Info

Fabric Choices



### 7408 Options: brushed aluminum base †

Series	Standard Height	with Crescent Back	with Backrest
7408-xx	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
7408CB-xx	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
7408B-xx	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>

★ options available at additional cost

### Specifications:

Series	Seat Diameter		Base Diameter with casters		Standard Seat Height				Overall Height <sup>w/</sup> Backrest				Caster Diameter hard wheel		Shipping Weight <sup>§</sup>		Weight Capacity <sup>†</sup> EVENLY DISTRIBUTED	
	inch	cm	inch	cm	low		high		low		high		inch	cm	lb	kg <sup>†</sup>	lb	kg <sup>†</sup>
7408	16.00	40.64	24.00	60.96	19.50	49.53	26.50	67.31	35.50	90.17	42.50	107.95	2.00	5.08	45	20.41	350	158.76

★ options available at additional cost

† -xx in the product number is the placeholder for fabric choice

† standard rounding method used

<sup>†</sup> maximum weight capacity - evenly distributed - refer to user manual

<sup>§</sup> approximate measurement



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF24100055 Rev A24