

STETHOSCOPES


LIGHTWEIGHT SINGLE HEAD STETHOSCOPE

The Lightweight Single Head Stethoscope is ideal for general assessment and blood pressure measurement. The slim, lightweight anodized aluminum chestpiece easily slides beneath a blood pressure cuff. Graham-Field stethoscope headsets are anatomically adjustable to ensure a comfortable acoustically sealed aural fit.

Features

- Complete unit weighs only 2.5 oz (70 g)
- 1 3/4" (4.5 cm) diameter diaphragm with non-chill diaphragm ring
- Adjustable anodized aluminum binaurals
- Flexible PVC tubing
- 22" (56.5 cm) tubing; 31.5" (80 cm) total length
- Full color retail packaging
- Includes stethoscope with diaphragm only chestpiece, adjustable binaurals, and (1) set of white firm threaded PVC eartips
- This product is not made with natural rubber latex
- One (1) year limited warranty

Item	Color		UOM
300DLX	Black		1 ea
300DLX-BL	Blue		1 ea
300DLX-GY	Grey		1 ea
300DLX-R	Red		1 ea
300DLX-Y	Yellow		1 ea
300DLX-LB	Light Blue		1 ea
300DLX-TL	Teal		1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.