

## STETHOSCOPES

# CARDIOLOGY STETHOSCOPE

### Features

- Features a stainless steel combination chestpiece for maximum convenience and internal spring binaurals for high acoustic sensitivity
- Dual non-chill rings provide patient comfort and 22" (56.5 cm) latex-free tubing keeps acoustics clear, even in the noisiest settings
- Comes complete with additional diaphragm and two pair of contoured eartips



Item	Color	UOM
425	Grey	1 ea
425BLK	Black	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.