



## HEALTHCARE SEATING



# PHLEBOTOMY CHAIR

### 6200 SERIES

The Lumex 6200 Series bariatric blood draw/phlebotomy chair is designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

#### Features:

- Fully upholstered padded seat, back, and arm rests
- Shown with optional wall saver wireframe
- Powder coat frame finish - Foghorn
- Available in 24 industry standard colors to complement most decors
- Individual height adjustable arm rests

[6200 Product Info](#)

[Fabric Choices](#)



Single Flip Arm  
Dual Flip Arm  
Wall Saver Wireframe

#### 6200 Options:

| Series                | Single Flip Arm | Dual Flip Arm | Wall Saver Wireframe |
|-----------------------|-----------------|---------------|----------------------|
| 6200-xx               | ■               |               |                      |
| 6200-DA-xx            |                 | ■             |                      |
| 6200-WALL SAVER-xx    | ■               |               | ■                    |
| 6200-DA-WALL SAVER-xx |                 | ■             | ■                    |

★ options available at additional cost



#### Specifications:

| Series | Overall Width |       | Overall Depth |       | Overall Height |       | Seat Width |       | Seat Depth |       | Seat Height |       | Armrest Height |       | Back Width |       | Back Height |       | Shipping Weight <sup>2</sup> |                 | Weight Capacity <sup>1</sup><br>EVENLY DISTRIBUTED |                 |     |        |
|--------|---------------|-------|---------------|-------|----------------|-------|------------|-------|------------|-------|-------------|-------|----------------|-------|------------|-------|-------------|-------|------------------------------|-----------------|--|-----------------|-----|--------|
|        | inch          | cm    | inch          | cm    | inch           | cm    | inch       | cm    | inch       | cm    | inch        | cm    | low            | high  | inch       | cm    | inch        | cm    | lb                           | kg <sup>3</sup> | lb   | kg <sup>3</sup> |     |        |
| 6200   | 34.00         | 86.36 | 29.00         | 73.66 | 37.50          | 95.25 | 24.50      | 62.23 | 18.50      | 46.99 | 22.00       | 64.77 | 27.00          | 62.23 | 36.00      | 66.04 | 25.50       | 64.77 | 17.00                        | 43.18           | 55   | 24.95           | 800 | 362.87 |

★ options available at additional cost † -xx in the product number is the placeholder for fabric choice † standard rounding method used  
<sup>1</sup> maximum weight capacity - evenly distributed - refer to user manual <sup>2</sup> overall height without backrest <sup>3</sup> approximate measurement



GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
 The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
 Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF1800043 Rev G24