


STETHOSCOPES

PEDIATRIC STETHOSCOPE

Features

- Lightweight anodized aluminum rotating dual head chestpiece with recessed non-chill ring
- Adjustable chrome-plated brass external spring binaural
- 22" (56.5 cm) latex-free tubing



Item	Color	UOM
507	Black	1 ea
507LB	Light Blue	1 ea
507P	Pink	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.