

STETHOSCOPES

TEACHING STETHOSCOPE

Features

- Excellent acoustical clarity
- 22" (56.5 cm) latex-free tubing with connector allows both student and instructor to listen simultaneously
- Anodized aluminum combination dual head chestpiece and chrome-plated brass rotating stem
- Adjustable chrome-plated external spring binaural



Item	Color	UOM
540	Black	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.