

STETHOSCOPES

SPRAGUE RAPPAPORT-TYPE 22" STETHOSCOPE

Features

- Five interchangeable chestpieces (two diaphragms, three bells) facilitates detection of high, medium, and low frequency sounds
- Ideal for infant, pediatric and adult use
- Includes three ear tip sets to maximize comfort
- Steel binaural
- Heavy-walled 22" (56.5 cm) tubing decreases extraneous sounds
- Accessory pouch included
- Latex free



Item	Color		UOM
602	Black		1 ea
602BL	Blue		1 ea
602R	Red		1 ea
602TL	Teal		1 ea
602LB	Light Blue		1 ea
602MAG	Magenta		1 ea
602L	Lavender		1 ea
602BUR	Burgundy		1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.