

TUBING AND CLAMPS

# SIDE OPENJAW HOFFMAN SCREW COMPRESSOR

**Features**

- Serrated thumb screw head
- Chrome plated copper



		overall Width		overall Height		for Tube Size		UOM
		inch	cm	inch	cm	inch	cm	
3083DZ	Standard	7/8	2.2	1	2.5	0.5	1.3	12/bx



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2023, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.