





# AUTOMATIC OVERBED TABLES




*Total  
Design  
Solution*



### FINISHES AVAILABLE ON 3DL THERMOFORM OVERBED TABLE TOPS




Aria



Cocobala



Kensington Maple



Shaker Cherry




Solar Oak



Talas Cherry\*\*  
*\*\*Premium upcharge may apply*




Wild Cherry




Windsor Mahogany


Color variations may occur due to monitor or printer settings.



### FINISHES AVAILABLE ON LAMINATE OVERBED TABLE TOPS




Cocobala




Aria



Fonhill Pear



Kensington Maple



Shaker Cherry




Solar Oak



Talas Cherry  
*\*\*Premium upcharge may apply*



Wild Cherry



Windsor Mahogany

Color variations may occur due to monitor or printer settings.

AUTOMATIC OVERBED TABLES

Ready Ship Sets: Overbed Table Tops & Base <sup>5</sup>

	Item No.	A793053	A793020	A793029	A793004	A860053	A860029	A860004	A797053	A797020	A797029	A797004	A859053	A859029	A859004	Oversized Tops			
		A919053XL	A919020XL	A919029XL	A919004XL														
Low Height <sup>1</sup>	base height	■	■	■	■	■	■	■	—	—	—	—	—	—	—	—	—	—	—
Standard Height <sup>1</sup>	base height	—	—	—	—	—	—	—	■	■	■	■	■	■	■	—	—	—	—
Extended Height <sup>1</sup>	base height	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—
H Base	base style	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—
Modified H Base	base style	■	■	■	■	■	■	■	■	■	■	■	■	■	■	—	—	—	—
material <sup>3DL</sup> <sup>2</sup>		■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
Flat Edge	top profile	■	■	■	■	—	—	—	■	■	■	■	—	—	—	—	—	—	—
Raised Edge	top profile	—	—	—	—	■	■	■	—	—	—	—	■	■	■	—	—	—	—
Solar Oak	top finish	■	—	—	—	■	—	—	■	—	—	—	■	—	—	—	—	—	—
Cocobala	top finish	—	■	—	—	—	—	—	—	■	—	—	—	—	—	—	—	—	—
Windsor Mahogany	top finish	—	—	■	—	—	■	—	—	—	■	—	—	■	—	—	—	—	—
Wild Cherry	top finish	—	—	—	■	—	—	■	—	—	—	■	—	—	—	—	—	—	■
Weight Capacity EVENLY DISTRIBUTED	lb	40	40	40	40	40	40	40	40	40	40	40	40	40	40	100	100	100	100
UOM	ea	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1

■ Standard Feature ■ Optional Feature ■ Available Style ■ Available Finish — Not Available/Does Not Apply

Made To Order:

Top Options <sup>7</sup>

Item	Top Description	3DL <sup>2</sup> material	Laminate <sup>3</sup> inch	Width inch	Length inch	Raised Edge top profile	Flat Edge	Finish	UOM
A28B	Laminate Top	—	■	15	30	—	■	Specify <sup>3,4</sup>	1/ea
A28-3DLS	3DL OB Top	■	—	15	30	—	■	Specify <sup>2,4</sup>	1/ea
A28-RES	3DL Raised Edge	■	—	15	30	■	—	Specify <sup>2,4</sup>	1/ea
A28A	Lam. Oversized	—	■	18	35	—	■	Specify <sup>3,4</sup>	1/ea
OS3DLOBTR	3DL Oversized <sup>8</sup>	■	—	18	35	■	—	Specify <sup>3,4</sup>	1/ea
OS3DLOBTF	3DL Oversized <sup>8</sup>	■	—	18	35	—	■	Specify <sup>3,4</sup>	1/ea

Base Options <sup>7</sup>

Item	Height	Width inch	Length inch	Height Range <sup>1</sup>		Weight Capacity EVENLY DISTRIBUTED lb	UOM
				inch - low	inch - hi		
A793000	Low	15.30	28.50	21.00	32.00	40	1/ea
A797000	Standard	15.30	28.50	27.75	40.00	40	1/ea
A919000	Extended	18.20	29.30	29.90	43.00	100	1/ea

Accessories <sup>6</sup>

Item	Description	UOM
A28VAN	Vanity with Pop Up Mirror [FITS ALL]	1/ea

<sup>1</sup> height adjusts infinitely  
<sup>2</sup> 3DL finish swatches on page 2  
<sup>3</sup> laminate finish swatches on page 3  
<sup>4</sup> must specify OB top finish when ordering

<sup>5</sup> Ready Ship tops & bases sold together  
<sup>6</sup> available for either Ready Ship or Made to Order  
<sup>7</sup> Made to Order tops and bases each sold separately  
<sup>8</sup> Requires A919000 Extended Height 100 lb base

• additional styles of OB top may be available upon request  
 • all measurements include a manufacturing variance of ± 0.25 inch  
 • cover images - left: A859004 40 lb OB top right: A919055 100 lb OB top



GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
 The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
 Basic American Medical Products is a trademark of GF Health Products, Inc.



GF120004RevK24