

# AUTOMATIC OVERBED TABLE

## 100 LB AUTOMATIC OBT BASE \*

### Features

- Supports up to 100 lb — EVENLY DISTRIBUTED<sup>1</sup>  
*\*100 lb base only compatible with oversized 18" x 35" top*
- Automatic/infinite stops within height travel range
- Modified "H" base allows use with chair or bed
- Low profile design allows overbed table to be used with limited clearance beds
- Tops are customizable in different finishes and edge options

### Made To Order: Top Finish Options



Item No.	Height	Width		Length		Height Range <sup>†</sup>		Weight Capacity <sup>1</sup> EVENLY DISTRIBUTED	UOM
		inch	cm	inch	cm	inch - low	inch - hi		
A919000	Extended	18.2	29.3	29.9	43.0	100	1 ea		



### Made To Order: Top Edge Options

Item No.	Construction	Width		Length		Raised Edge	Flat Edge
		inch	cm	inch	cm	top profile	top profile
OS3DLOBTR	3DL Oversized	18	35	—	—	■	—
OS3DLOBTF	3DL Oversized	18	35	—	—	—	■

### Ready Ship Sets: Oversized Overbed Table Tops & Base

Item No.	Oversized Top Finish	Oversized Top				100 lb Base				OBT Height Range				Weight <sup>1</sup> Capacity EVENLY DISTRIBUTED		
		Width		Length		Width		Length		minimum		maximum		lb	kg	UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm			
A919053XL	Solar Oak	18.00	45.72	35.00	88.90	18.20	46.23	29.30	74.42	29.90	75.95	43.00	109.22	100	45.36	1 set
A919020XL	Cocobala	18.00	45.72	35.00	88.90	18.20	46.23	29.30	74.42	29.90	75.95	43.00	109.22	100	45.36	1 set
A919029XL	Windsor Mahogany	18.00	45.72	35.00	88.90	18.20	46.23	29.30	74.42	29.90	75.95	43.00	109.22	100	45.36	1 set
A919004XL	Wild Cherry	18.00	45.72	35.00	88.90	18.20	46.23	29.30	74.42	29.90	75.95	43.00	109.22	100	45.36	1 set

GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.grahamfield.com](http://www.grahamfield.com)

