

## SPHYGMOMANOMETERS

# LABSTAR™ DELUXE SPHYGMOMANOMETER

### Features

- Accumax™ nylon cuff with index & range markings
- Latex two (2) tube bladder
- Large knurled air release valve for precision deflation
- Deluxe end valve reduces dust build-up
- Color coordinated gauge and cuff in blue, black, or grey
- Pressure measurement from 20 mmHg to 300 mmHg
- Lifetime calibration warranty
- Zippered carrying bag is included

Caution: This product contains natural rubber latex, which may cause allergic reactions.



Item	Size	Tubes	Bladder	Color	UOM
202	Adult	2-Tube	Latex	Blue	1 ea
202C	Child	2-Tube	Latex	Blue	1 ea
202I	Infant	2-Tube	Latex	Blue	1 ea
202N	Newborn	2-Tube	Latex	Blue	1 ea
202T	Thigh	2-Tube	Latex	Blue	1 ea
202X	Large Adult	2-Tube	Latex	Blue	1 ea
202BK	Adult	2-Tube	Latex	Black	1 ea
202GY	Adult	2-Tube	Latex	Grey	1 ea

	Limb Range		cuff Width		cuff Height	
	inch	cm	inch	cm	inch	cm
Newborn	3.75 – 5.25	9.5 – 13.3	7	17.8	2	5.1
Infant	5 – 7.5	12.7 – 19.1	10.5	26.7	3	7.6
Child	7 – 10.5	17.8 – 26.7	13.5	34.3	4	10.2
Adult	9.5 – 15	24.1 – 38.1	20.5	52.1	5.5	14
LG Adult	13 – 19	33 – 48.3	25	63.5	6.5	16.5
Thigh	17 – 26	43.2 – 66	29.5	74.9	8.5	21.6



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2011, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.