

SPHYGMOMANOMETERS

WALLMAX™ ANEROID SPHYGMOMANOMETER

Features

- Luminescent dial
- Large, 6" (15.2 cm), easy-to-read gauge
- Wall mount bracket with 110° swivel range
- Accumax™ nylon cuff with 8' (2.4 m) coiled tubing
- Cuff can be stored behind gauge
- One (1) year warranty for manufacturing defects only

Caution: This product contains natural rubber latex, which may cause allergic reactions.



Item	Size	Bladder	Color	UOM
222	Adult	Latex	Blue	1 ea
V222B	Adult	Non-Latex	Black	1 ea

	Limb Range		cuff Width		cuff Height	
	inch	cm	inch	cm	inch	cm
Newborn	3.75 – 5.25	9.5 – 13.3	7	17.8	2	5.1
Infant	5 – 7.5	12.7 – 19.1	10.5	26.7	3	7.6
Child	7 – 10.5	17.8 – 26.7	13.5	34.3	4	10.2
Adult	9.5 – 15	24.1 – 38.1	20.5	52.1	5.5	14
LG Adult	13 – 19	33 – 48.3	25	63.5	6.5	16.5
Thigh	17 – 26	43.2 – 66	29.5	74.9	8.5	21.6



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.