

SPHYGMOMANOMETERS

ACCUMAX™ COMPLETE INFLATION SYSTEMS

Features

- Calibrated nylon Accumax™ inflation cuff
- Two-tube bladder
- Deluxe bulb & valve



Item	Size	Tubes	Bladder	Color	UOM
2361A	Adult	2-Tube	Latex	Blue	1 ea
2361N	Newborn	2-Tube	Latex	Blue	1 ea
2361C	Child	2-Tube	Latex	Blue	1 ea
2361X	Large Adult	2-Tube	Latex	Blue	1 ea
2361T	Thigh	2-Tube	Latex	Blue	1 ea
2365A	Adult	2-Tube	Latex	Black	1 ea
2365X	Large Adult	2-Tube	Latex	Black	1 ea
2365T	Thigh	2-Tube	Latex	Black	1 ea
2360A	Adult	2-Tube	Non-Latex	Black	1 ea
2360C	Child	2-Tube	Non-Latex	Black	1 ea
2360T	Thigh	2-Tube	Non-Latex	Black	1 ea
2360X	Large Adult	2-Tube	Non-Latex	Black	1 ea
2366A	Adult	2-Tube	Non-Latex	Grey	1 ea
2367A	Adult	2-Tube	Non-Latex	Blue	1 ea
2367C	Child	2-Tube	Non-Latex	Blue	1 ea

	Limb Range		cuff Width		cuff Height	
	inch	cm	inch	cm	inch	cm
Newborn	3.75 – 5.25	9.5 – 13.3	7	17.8	2	5.1
Infant	5 – 7.5	12.7 – 19.1	10.5	26.7	3	7.6
Child	7 – 10.5	17.8 – 26.7	13.5	34.3	4	10.2
Adult	9.5 – 15	24.1 – 38.1	20.5	52.1	5.5	14
LG Adult	13 – 19	33 – 48.3	25	63.5	6.5	16.5
Thigh	17 – 26	43.2 – 66	29.5	74.9	8.5	21.6



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.