



SUPPORT SURFACES



FOAM MATTRESS

BARIATRIC CARE SERIES HIGH-DENSITY FOAM MATTRESS

Features

- The Bariatric Care Foam Mattress offers comfort, support and durability for the bariatric patient.

Cover Features

- Ultra soft nylon cover. Durable Vinyl cover bottom.
- Fluid resistant
- Anti-static and flame retardant
- Vapour permeable (0% vapor transmission)
- Infection Control is aided by Antimicrobial properties; bacteria resistant, anti-fungal, and non-allergenic
- Concealed full zipper
- Easy to clean and maintain

Tri-Laminate³ Mattress Construction

- High-density foam supportive base
- Secondary layer provides support and comfort
- Top layer is designed with die-cut foam to provide comfort and allows air to flow and cool the sleep surface
- Firm side perimeter construction

Credentials

- Designed to meet the requirements of HCPCS Code E0184
- Meets 16CFR1633, Made in USA

Warranty

- 60 month Mattress
- 12 month Cover

Product Info



DEFINED PERIMETER
aids in fall prevention during ingress / egress



SOFT HEEL ZONE
for patients with sensitive feet

Shipped compress rolled for cost effective freight

Largest Box Size:²
14 x 14 x 46 inches



Specifications

Item No.	overall ¹ Width		overall ¹ Length		overall Thickness		Product ⁸ Weight		Weight ¹ Capacity EVENLY DISTRIBUTED		Tri-Laminate ³ die-cut foam top layer Dual Layer soft high-density topper	Raised Side Bolsters	Defined Perimeter	Soft Heel Zone	Soft Heel Slope	Pressure Reduction	Zippered Cover	Ships Compressed/Rolled ²	
	inch	cm ¹	inch	cm ¹	inch	cm ¹	lb	kg ¹	lb	kg ¹									
GEN3980A	39	99.06	80	203.20	7	17.78	33	14.97	800	362.87	■	—	—	■	■	—	■	■	■
GEN3980B	35	88.90	80	203.20			35	15.88	1000	453.59	■	—	—	■	■	—	■	■	■
GEN4280A	42	106.68	80	203.20	7	17.78	35	15.88	800	362.87	■	—	—	■	■	—	■	■	■
GEN4280B	37	94.00	80	203.20			37	16.78	1000	453.59	■	—	—	■	■	—	■	■	■
GEN4880A	48	121.92	80	203.20	7	17.78	40	18.14	800	362.87	■	—	—	■	■	—	■	■	—
GEN4880B	42	106.68	80	203.20			42	19.05	1000	453.59	■	—	—	■	■	—	■	■	■
GEN5480A	54	137.16	80	203.20	7	17.78	45	20.41	800	362.87	■	—	—	■	■	—	■	■	—
GEN5480B	47	119.41	80	203.20			47	21.32	1000	453.59	■	—	—	■	■	—	■	■	■
GEN4884H	48	121.92	84	213.36	7	17.78	36	16.33	650	294.84	—	—	—	—	—	—	—	—	—
650-1501-1633N	42	106.68	80	203.20			27	12.25	600	272.16	—	—	—	—	—	—	—	—	—

¹ foam tolerance ±1 inch ² standard rounding method used ³ maximum weight capacity - evenly distributed - refer to user manual ⁴ approximate measurement ⁵ box length depends on mattress width: 40" box for 36" mattress and 46" for 42". ⁶ item numbers ending with A or B



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.