



SUPPORT SURFACES



CLASSIC FOAM MATTRESS

31 SERIES DUAL LAYER FOAM MATTRESS

Features

- Provides an excellent value for budget-conscious facilities
- Specifically developed with patient comfort and pressure injury management in mind
- Provides pressure reduction and redistribution of stage 1 pressure ulcers

Cover Features

- Ultra soft nylon cover. Durable Vinyl cover bottom
- Vapour permeable (0% vapor transmission)
- Fluid resistant
- Infection control is aided by antimicrobial properties; Bacteria resistant, anti-fungal, and non-allergenic
- Anti-static and flame retardant
- Concealed full zipper (319 series only)
- Easy to clean and maintain

Mattress Construction

- Dual layer foam construction with convoluted top layer
- Convoluted top layer provides comfort and allows air to flow and cool the sleep surface



Credentials

- Designed to meet the requirements of HCPCS Code E0184
- Meets 16CFR1633, Made in USA

Warranty

- 24 month Mattress
- 12 month Cover

Product Info



CONVOLUTED TOPPER creates excellent airflow

Shipped compress rolled for cost effective freight

Largest Box Size:²
14 x 14 x 46 inches



Specifications

Item No.	Overall [†] Width		Overall [†] Length		Overall Thickness		Product Weight [§]		Weight ¹ Capacity EVENLY DISTRIBUTED		Convoluted Topper	Raised Side Bolsters	Defined Perimeter	Soft Heel Zone	Soft Heel Slope	Pressure Reduction	Zippered Cover
	inch	cm [†]	inch	cm [†]	inch	cm [†]	lb	kg [†]	lb	kg [†]							
31680-1633			80	203.20			17	7.71			■	—	—	—	—	■	■
31980-1633							17	7.71			■	—	—	—	—	■	■
31684-1633	36	91.44	84	213.36	6	15.24	18	8.16	300	136.08	■	—	—	—	—	■	■
31984-1633							18	8.16			■	—	—	—	—	■	■
31675-1633			75	190.50			16	7.26			■	—	—	—	—	■	■
31975-1633							16	7.26			■	—	—	—	—	■	■

[†] foam tolerance ±1 inch ¹ standard rounding method used ¹ maximum weight capacity - evenly distributed - refer to user manual
² box length depends on mattress width: 40" box for 36" mattress and 46" for 42". [§] approximate measurement



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF070011 Rev F24