



CLINICAL CARE RECLINERS

FR577RG SERIES MULTI-CARE RECLINER

Three reclining positions for patient comfort plus a caregiver-activated Trendelenburg (or Shock) Position.

Softer seat and back upholstery for improved patient comfort. Varied density seat back and redesigned seating position provide excellent comfort and support during long term sitting.

Improved wider and softer armrest for patient comfort.

Top section of headrest can be moved forward for enhanced patient comfort.

Easy to activate gas cylinder mechanism for quick Trendelenburg activation and improved lift.

Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Side table surface is smooth and features tapered edges for easy cleaning.

Easily removed side panels enable easy access to recliner for cleaning and maintenance.

Four Tente[®] swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster.

Wide selection of fabrics and colors. See the [Swatch Card](#) for details.



Articulating Headrest
in Forward Position

Options:

- Heat and Massage
- Heat Only feature
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow



Multiple Positions for Enhanced Patient Comfort

Model (xxxx represents fabric colors)	Wall Saver Armrest Edges	Fully Upholstered Arms	Heat & Massage	Heat Only	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Width between side tables in upright position		Depth fully reclined (Trendelenburg)		Weight Capacity EVENLY DISTRIBUTED		UOM
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
FR577RG	●				31.5	80	39	99.1	48	121.9	21.5	54.6	19.5	49.5	20.5	52.1	6.5	16.5	21.75	55.2	52	132.1	75	190.5	290	131.5	1 ea
FR577RGH	●		●		31.5	80	39	99.1	48	121.9	21.5	54.6	19.5	49.5	20.5	52.1	6.5	16.5	21.75	55.2	52	132.1	75	190.5	290	131.5	1 ea
FR577RGHO	●			●	31.5	80	39	99.1	48	121.9	21.5	54.6	19.5	49.5	20.5	52.1	6.5	16.5	21.75	55.2	52	132.1	75	190.5	290	131.5	1 ea
FR577RGxxxxPA		●			31.5	80	39	99.1	48	121.9	21.5	54.6	19.5	49.5	20.5	52.1	6.5	16.5	21.75	55.2	52	132.1	75	190.5	290	131.5	1 ea
FR577RGHxxxxPA		●	●		31.5	80	39	99.1	48	121.9	21.5	54.6	19.5	49.5	20.5	52.1	6.5	16.5	21.75	55.2	52	132.1	75	190.5	290	131.5	1 ea
FR577RGHOxxxxPA		●		●	31.5	80	39	99.1	48	121.9	21.5	54.6	19.5	49.5	20.5	52.1	6.5	16.5	21.75	55.2	52	132.1	75	190.5	290	131.5	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.

The most current and complete product information can be found on our website.

www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.