



EXTENDED CARE BEDS



# ZENITH SERIES

**7200 & 9200 + SLIDE-W-I-D-E®**



The Zenith™ Series by Basic American is the Premier Extended Care Bed. Engineered for durability and function, The Zenith features sealed ball bearings in the joints for smooth, quiet movement and SilverSolutions in the powder coating to provide antimicrobial surface protection.

Basic American by Graham-Field is a leading manufacturer of beds serving the senior living industry for 40+ years. We proudly manufacture our beds in the USA, in Fond du Lac, Wisconsin, which enables faster response to customers' needs.

Innovation with industry-leading technology is designed into the Basic American Zenith, Matrix, and Liberty models to provide feature-rich beds and value for every budget.

**7200**  
ZENITH APS BED  
*shown with pivot assist bar (right),  
counter-rotating assist (left) device  
available to add anytime at additional cost*

## ZENITH SERIES

### FEATURES AND ADVANTAGES:

#### Made in the USA

- Minimal lead times
- Quality craftsmanship and manufacturing

#### SilverSolutions

- Located in the powder-coating which provides antimicrobial surface protection
- Premium Products

#### Sealed Ball Bearings

- Allows seamless, continuous, and noiseless movement
- Maintenance Free

#### SLIDE-W-I-D-E® Integrated Expandable Deck

- Slides located on the frame effortlessly slide in and out allowing for quick and easy transition
- Expandable/retractable deck adjusts to 35", 39", and 42" with two quick-release ball-detent pins on each side

#### ANSI/AAMI ES60601 Safety Certification

- Tip stability testing
- Load capacity testing (twice the operating capacity)
- Roll clearance testing
- Splash proof electronics testing

#### Optional 'Plug and Play' digital electronics enable future feature additions

- Staff Control for Trendelenburg and Reverse Trendelenburg
- Chair Position
- Underbed Light
- USB Charger

#### Certified as a transport device for the CMS Emergency Preparedness Plan

- Easy maneuvering/adjustability through doorways and hallways during emergency or evacuation

#### Intuitive two-pedal locking system

- Centrally located
- Easily visible and accessible at ALL heights and widths

#### Industry-leading, energy-saving electronics

- Reduces Energy Consumption

#### Grounded bed with right-angle hospital grade plug

- Provides added safety for both staff and residents

#### US Patented design of high/low mechanism


- Housed with the mainframe tubes; no holes are located on the top side that will harbor bacteria

#### Grid sleep surface with cross braces

- Sturdy, no-sag sleeping surface

### Warranty:

Main Frame & Welds	15 years
Control Box & Actuator Motors	5 years
Hand Control Pendant & Cabling	3 years
Staff Control & Cabling	3 years
Head & Footboards	1 year



**7200**  
ZENITH™ APS BED WITH  
**Slide-W-I-D-E**  
EXPANDABLE/RETRACTABLE SLEEP DECK



USB power supply (right)  
underbed light (below)



**9200**  
ZENITH BED  
shown with pivot assist bar (right),  
counter-rotating assist (left) device  
available to add anytime at additional cost



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Basic American Medical Products and Slide-W-I-D-E are registered trademarks of GF Health Products, Inc.

graham  
field

GF2100055 RevE24

± 1 inch tolerance